

AV-S Ventilator

The intuitive user interface and comprehensive support modes provide the optimum therapy for all patient profiles



- ◆ Combines sophistication and ease of use
- ◆ Large, colour Touchscreen and Com-wheel control
- ◆ Volume and Pressure ventilation, plus SMMV, SIMV, PSV and PEEP
- ◆ Single/dual waveform display
- ◆ High quality, multi-option product with flexible specification
- ◆ 'Life-Care' or optional 'Life-Care Plus' Warranty and Customer Care scheme

Partnership for Life

AV-S Ventilator

Multi-mode Ventilator

An easy to use, multifunction anaesthesia ventilator, designed for all patient profiles

- ◆ Volume, PCV, PSV, SIMV and SMMV modes
- ◆ Comprehensive printer/data outputs for networking and interfacing to patient monitors
- ◆ Integrated Oxygen Monitor and Spirometry
- ◆ Inverse I:E Ratio capability
- ◆ Electronic PEEP
- ◆ Autoclavable Latex Free Bellows
- ◆ Oxygen or Air drive gas
- ◆ 30 minutes battery backup
- ◆ Selectable Dual Waveform Display:
 - Pressure v. Time
 - Volume v. Time
 - Pressure v. Volume (for ventilation analysis) plus waveform freeze facility
- ◆ Save and recall function for user specific settings
- ◆ Adult and Paediatric default settings
- ◆ Flexible specification
 - Stand-alone operation or seamless integration with Prima SP Workstation
 - Display mounting option
 - Multilingual display



AV-S Ventilator

Partnership for Life

Penlon's philosophy embraces commitment to a successful, long term relationship with all our customers




Touchscreen Control

Prompt and easy operation from either the Touchscreen or via the Com-wheel




Latex-free bellows

The latex-free bellows, canister and base assembly are fully autoclavable



SIMV

Combines spontaneous and set mandatory breaths



SMMV

Combines spontaneous breaths with synchronised mandatory breaths to achieve the set minute volume



PSV

Assists each spontaneous breath with a preset pressure, thus reducing the effort required to breathe

Advanced Spontaneous Breathing Modes


The AV-S ventilator provides three support modes that can be utilised as the patient attempts to breathe spontaneously. Patient recovery is accelerated by increased tidal volume and SpO₂, and reduced EtCO₂.

SIMV Synchronised Intermittent Mandatory Ventilation

Guarantees a minimum level of volume. SIMV allows spontaneous breaths and set mandatory breaths


SIMV - Spontaneously Breathing Patient

Negative pressure[†] in the Trigger Window* (generated by the patient's spontaneous breath) results in a synchronised mandatory breath at a preset volume and rate



SIMV - Non-breathing Patient

If the patient makes no effort to breathe during a cycle, a mandatory breath, at the end of the Trigger Window,* will still be delivered at the preset volume and rate




SMMV Synchronised Mandatory Minute Ventilation

Guarantees a set level of minute volume ventilation. SMMV allows spontaneous breaths, combined with a synchronised mandatory breath, to achieve the set minute volume


SMMV - Spontaneously Breathing Patient

Negative pressure[†] in the Trigger Window* (generated by the patient's spontaneous breath) results in a synchronised mandatory breath, ensuring that the set minute volume is achieved



SMMV - Non-Breathing Patient


If the patient makes no effort to breathe during a cycle, a mandatory breath, at the end of the Trigger Window,* will still be delivered at the preset minute volume and rate



PSV Pressure Supported Ventilation

PSV assists each spontaneous breath with a preset pressure, thus reducing the effort required to breathe. Negative pressure[†] (generated by the patient's spontaneous breath) results in synchronised pressure support

PSV is used to support spontaneously breathing patients ONLY If the patient makes no attempt to breathe, the ventilator will not provide support and the apnoea alarm will be activated



*Trigger Window = 60% of Cycle Time • [†]Negative relative to PEEP

Technical Specification

Physical	
Size	- Control Unit Only - with Adult Bellows
Screen	
Weight	- Control Unit Only - with Adult Bellows
Bellows (Latex free)	
Power	
Drive Gas	
Functional	
Tidal Volume (Vt)	
Rate (BPM)	
I:E Ratio	
Pressure Limit	
Fresh Gas Compensation	
Ventilation Modes	
Sigh Function (Volume Mode)	
Pressure Control	
Spontaneous Mode	
Electronic PEEP	
Oxygen Monitor	
SIMV, SMMV, PSV	
Trigger	
Trigger Window	
Tidal Volume (Vt)	
Minute Volume (Vm)	
Inspiratory Time (Ti)	
Support Pressure	

For further information please contact your local Penlon representative

Information contained in this leaflet is correct at the date of publication. The policy of Penlon Limited is one of continued improvement to its products. Because of this policy Penlon Limited reserves the right to make any changes, which may affect the information in this leaflet without giving prior notice.

Alarms - Automatic		
Alarm Mute	30 Seconds	
Low Drive Gas Pressure	Less than 235 kPa (34 psi)	
High Continuous Airway Pressure	Above 30 cmH ₂ O at start of cycle	
Low Pressure	4 to 14 cmH ₂ O PEEP Referenced	
Low Tidal Volume	50% of Volume Set (Spirometry)	
Incorrect Rate or Ratio		
Mains Failure	30 Minutes Battery Backup	
Low Battery	5 Minutes Use	
Vent Inop	Internal or Battery Failure	
Apnoea	Flow Referenced	
Alarms - Optional User Set		
Tidal Volume	- Minimum	0 to 1600 ml
	- Maximum	20 to 1600 ml
Minute Volume	- Minimum	0 to 10 L
	- Maximum	0 to 30 L
Low and High O ₂ Concentration	18% to 105%	
High Airway Pressure	10 to 80 cmH ₂ O Adjustable	
Default Settings	Adult	Paediatric
VOLUME		
• Tidal Volume (Vt)	600 ml	150 ml
• Rate (BPM)	10	15
• I:E Ratio	1:2	1:2
• Pmax	38 cmH ₂ O	38 cmH ₂ O
PRESSURE		
• Tidal Volume (Vt)	600 ml	150 ml
• Rate (BPM)	10	15
• I:E Ratio	1:2	1:2
• P-Target	10 cmH ₂ O	10 cmH ₂ O
SIMV		
• Tidal Volume (Vt)	600 ml	200 ml
• Rate (BPM)	6	10
• Inspiratory Time	2 Seconds	1 Second
• Trigger	-1 cmH ₂ O	-1 cmH ₂ O
SMMV		
• Minute Volume (Vm)	3.6 L	2 L
• Rate (BPM)	6	10
• Inspiratory Time	2 Seconds	1 Second
• Trigger	-1 cmH ₂ O	-1 cmH ₂ O
PSV		
• Support Pressure	10 cmH ₂ O	10 cmH ₂ O
• Inspiratory Time	2 Seconds	1 Second

For further information please contact
your local Penlon representative

Information contained in this leaflet is correct at the date of publication.
The policy of Penlon Limited is one of continued improvement to its products.
Because of this policy Penlon Limited reserves the right to make any changes,
which may affect the information in this leaflet without giving prior notice.



Penlon Limited Abingdon Science Park Barton Lane Abingdon OX14 3PH UK

www.penlon.com

General
Tel: +44 (0) 1235 547000
Fax: +44 (0) 1235 547041
http://www.penlon.com

International Sales
Tel: +44 (0) 1235 547001
Fax: +44 (0) 1235 547021
E-Mail: export@penlon.co.uk

UK Sales
Tel: 01235 547036
Fax: 01235 547023
E-Mail: uksales@penlon.co.uk

Service
Tel: 01235 547060
Fax: 01235 547061
E-Mail: service@penlon.co.uk

Medical Gas Solutions
Tel: +44 (0) 1235 547038
Fax: +44 (0) 1235 547055
E-Mail: medgas@penlon.co.uk