

ELECTRODES & ACCESSORIES

Maximum result both in terms of stability and precision, also under the most difficult conditions. The line of electrodes is extremely complete and makes it possible to meet any monitoring requirement.

- For adults or for paediatric use • With solid gel or liquid gel • For normal monitoring or for stress tests • FOAM or TNT/TEXTILE or MICROPORE base • Round-shaped with different diameters, oval or rectangular.

This line of accessories offers complementary or similar articles that are in great demand. Also in the case of accessories, the continuous research and exchange of information with the final user have allowed us to manufacture high-quality and extremely reliable products.



CERACARTA

...for over
50 years
 CERACARTA
 Italy
 since...1960

MONITORING & RESTING ECG'S



NM 50 RFI

Ø50 mm.	Dimensions	Ø55 mm.
Foam	Backing	Foam
Solid	Gel	Liquid
Ag/AgCl	Sensor	Ag/AgCl

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

Packaging

25 pcs / pouch - 500 pcs / box - 6000 pcs / carton

The most typical sizes. The most demanded supports. A proved quality always useful.



NM 55 RLI



NM 42 DFI

45x42 mm.	Dimensions	42x36 mm.
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

Packaging

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

Two multifunctional models: suitable for adult and pediatric ECG's.



NM 3642 OFI



NP 30 RFI

Ø30 mm.	Dimensions	28x44 mm.
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

Packaging

30 pcs / pouch - 600 pcs / box - 7200 pcs / master carton

The classical and largely used electrode for pediatric ECG's.
Excellent quality ensured also on "critical" patients.



NM 2844 RFS



MEDTAB

23x34 mm.	Dimensions	28x44 mm.
Silicone coated pet film	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl - Banana plug

100 pcs / pouch - 500 pcs / box - 5000 pcs / carton

Packaging

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

This electrode entirely replaces suction chests and clamps to patient's advantage both in terms of hygiene and health.



NM 2844 FBI

STRESS TEST



ST 50 RFI

Ø 50 mm.	Dimensions	Ø 50 mm.
Foam	Backing	Foam
Solid	Gel	Liquid
Ag/AgCl	Sensor	Ag/AgCl



ST 50 RLI

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

Packaging

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

The most suitable electrodes for ECG's under difficult conditions.
Easy to apply and use • Strong adhesiveness • Do not irritate the skin • Highly conductive gel

HOLTER/LONG MONITORING



NM 55 TNI

Ø 55 mm.	Dimensions	Ø 50 mm.
TNT	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Carbon

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

Packaging

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

When a larger adhesive surface is mandatory for an excellent examination. Easy to apply and use • Strong adhesiveness • Do not irritate the skin • Highly conductive gel

RADIO-TRANSPARENT



RT 50 RFI

MONITORING / PRE-WIRED ELECTRODES

RE 30 LWF

30 mm.
female 4 mm

	Dimensions	
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

3 pcs / pouch - 300 pcs / box

Packaging

RE 22 LWS

22 mm.
safety 1,55 mm

3 pcs / pouch - 300 pcs / box

Available also with 30 mm diameter Model RE 30 LWS (with 1,55 mm safety plug)

RT 30 LWF

30 mm.
female 4 mm

	Dimensions	
Foam	Backing	Foam
Solid	Gel	Solid
Carbon	Sensor	Carbon

3 pcs / pouch - 300 pcs / box

Packaging

RT 22 LWS



22 mm.
safety 1,55 mm

3 pcs / pouch - 300 pcs / box

Available also with 30 mm diameter Model RT 30 LWS (with 1,55 mm safety plug)



Available also with 22 mm diameter Model RE 22 LWF (with 4 mm female plug)



Available also with 22 mm diameter Model RT 22 LWF (with 4 mm female plug)



MEDI SENSOR LINE

MONITORING & RESTING ECG'S

MODEL F 50-S

MODEL F 3642-S



Ø 50 mm.	Dimensions	Ø 42.5 x 36.5 mm.
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

30 pcs / pouch – 1200 pcs / carton

Available also with liquid gel Model F 50-L

Packaging

50 pcs / pouch – 2000 pcs / carton


Available also with liquid gel Model F 3642-L



HOLTER / LONG MONITORING

MODEL W 601-S

MODEL T 601-S



Ø 55 mm.	Dimensions	Ø 55 mm.
TNT	Backing	MICROPORE
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

30 pcs / pouch – 1500 pcs / carton

Available also with Ø 50 mm Model W 601
and also with liquid gel Model W 601-L

Packaging

30 pcs / pouch – 1500 pcs / carton



STRESS TEST

MODEL FS 50-L

MODEL F 55-L



Ø 50 mm.	Dimensions	Ø 55 mm.
Foam	Backing	Foam
Liquid	Gel	Liquid
Ag/AgCl	Sensor	Ag/AgCl

30 pcs / pouch – 1200 pcs / carton

Packaging

30 pcs / pouch – 1200 pcs / carton



MODEL FS 501-S

MODEL F 50-S BLUE




Ø 35x45 mm	Dimensions	Ø 50 mm.
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Carbon

50 pcs / pouch – 2000 pcs / carton

Packaging

30 pcs / pouch – 1200 pcs / carton

Ideal for all ECG monitoring, radiolucent application included!




SUCTION CHEST ELECTRODES



Ag/AgCl sensor to guarantee maximum reliability in signal transmission.

- Available in three sizes: Ø15 mm for pediatric use / Ø 24 mm and Ø 30 mm for adult use.

RUBBER LIMB & THORACIC STRAPS



For attaching electrodes • for limbs use - one hole [15 mm x 40 cm / 25 mm x 45 cm], - two holes [32mm x 45cm] - for thoracic use - one hole [32mm x 135cm], six holes [100mm x 135cm].

CLAMP ELECTRODES



Classical-designed clamps with an improved quality. Ag/AgCl to ensure excellent signal transmission • Silver component to prevent distortions and to increase efficiency • Universal connection.

ELECTRO SURGICAL PLATES



ELECTROSURGICAL PLATES WITH A SUPERIOR QUALITY

- Made of FOAM suitable to obtain perfect results with all instruments • manufactured and tested in compliance with EN 60601-1/EN 60601-2/AAMI-HF 18 • available with one zone and two zones and REM alarm for adults (adhesive surface 215 sq cm), for children (surface 129 sq cm) and for new born (surface 63 sq cm). **All the models are available even equipped with cables.**

DEFIBRILLATOR PADS



A wide range of pads suitable for all the defibrillators available. Raw material quality and attention to details resulted into an absolutely reliable range of products. Universal without cable or equipped with cable, having a specific connection to comply with all instruments.

ELECTRODES FOR TENS AND EMG




• Superior quality • available with single or double cable • lead wire with 2 mm socket or press-stud connection • manufactured in compliance with EN ISO 9001 and AAMI standards.

PULSE OXIMETER SENSORS



A broad range of sensors to monitor and detect the arterial oxygen saturation. Excellent quality and proven accuracy are the characteristics that better describe a whole series of disposable and reusable sensors. Suitable for all the main fittings spread on the market.

ELECTRODES & ACCESSORIES FOR EEG & EMG



A) BRIDGE ELECTRODES Made of silver/silver chloride to optimize signal transmission. Slightly curved support to improve head adaption.

B) SILICONE HEADCAPS With 2 or 4 adjustable bands to fit for various requirements. Anatomically-designed earmuffs. Available in three sizes [Small-Medium-Large];

C) SPECIAL INK FOR EEG Suitable & excellent for all devices [Siemens excluded], allowing a perfect trace without blurs. Black colour- 1 litre bottle

D) DISC ELECTRODES Available in silver/silver chloride or gold - 10 mm diameter- With 120-150-200-250 cm lead wire & DIN 42802 safety plug - 10 electrodes per set with coloured lead wires (available 10 colours);

E) DISPOSABLE CONCENTRIC NEEDLES FOR ELECTROMYOGRAPHY Ergonomic handgrip and extremely sharp point - Equipped with wings for automatic unfastening -Shielded, light & flexible lead wire with 5 pole DIN 45322 connector and entirely shielded handgrip. Different diameters (from 0.35 to 0.60 mm) and different lengths (from 20 to 70 mm) - 50 pcs per packing.

MOUTHPIECES



Made of material particularly resistant to saliva or manipulation, perfectly hygienic, without residues and non-allergic • can be used on most spirometers of the most important brands present on the market • available in the external diameters 22-26-28-30 and 34 mm, but in case of large orders the customer can request a specific diameter • boxes of 500 pieces • an extremely handy dispenser is available on customer's demand.



Via Secondo Casadei N.14
zona industriale Villa Selva - 47122 Forlì (Italy) tel. +39.0543.780055 - fax +39.0543.781404
www.ceracarta.it - e mail: export@ceracarta.it