



GELS & CREAMS

SUPERIOR QUALITY

Ceracarta produces Gels for ultrasounds, for ECG and for EEG whose quality can without any doubt be defined as superior. The results obtained in this field come from a detailed research, combined with several decades' experience as well as continuous tests and trials in laboratories and in the most modern hospitals. All this has convinced even the most demanding users that the reliability of Ceracarta Supergels is unbeatable and its products are regarded among the best found on the market today.



CERACARTA





ECO SUPERGEL FOR ULTRASOUNDS

ECO SUPERGEL FOR ULTRASOUNDS

The carefully designed composition of this gel means that it is particularly suitable for diagnostic exams and ultrasound treatments.

- No salt is contained to prevent damaging the probes
- It does not cause skin irritation even after many scans
- Excellent viscosity (80,000 RPM5-STD 18")
- Rapid transmission speed (1.48)
- Perfect acoustic impedance
- Extremely efficient
- Bacteriostatic levels lower than admissible even by the most stringent international standards

AVAILABLE IN THREE VERSIONS:

- 260 g (box with 25 bottles)
- 1 kg (box with 8 bottles)
- 5 kg (individual pack with bottle/refill - box with 4 pieces)



Available in
"CLEAR"
version too.



ECG SUPERGEL & SUPERCREAM



ECG SUPERGEL & SUPERCREAM

In order to obtain high-quality ECG recordings for a long time it is essential to have a really "super" gel. The composition studied for the gel or for the cream (same characteristics, but with a lightly perfumed cosmetic base) actually allows for an improvement of the signal, preventing any kind of anomaly and providing users with a constantly reliable monitoring.

- No salt contained
- It does not damage the electrodes it comes into contact with
- Soluble in water
- It does not irritate even the most delicate skin types
- Not greasy
- Optimal viscosity
- Bacteriostatic levels lower than admissible even by the most stringent international standards
- PP-HDPE 260 g container with total guarantee (box with 25 bottles)

SAFE-E-SPRAY



SAFE-E-SPRAY / ECG & EEG

An absolutely practical spray to obtain the best possible conductivity during ECG/EEG monitoring. It strengthens the performance of the electrodes, by removing the bacteria and the dirt gathered on the surface. This product is recommended for preparing the skin before examination. Hypoallergenic, water soluble, antibacterial, odourless, non-greasy, does not stain. Packing: 250 ml. bottle ; n°25 bottles/box.

SUPER CLEAN



SUPER CLEAN

Degreasing paste apt for skin cleansing before EEG-EP-ECG recordings.

Allowing to reduce impedance down to minimum levels when a signal devoid of any adulteration is needed.

- Water-based, viscous gel
- Colouring agents free, non irritating
- 160 g. PLTHD-PLTLD bottle (24 pcs per box)



EEG/EMG SUPERCREAM & PASTE

EEG SUPERCREAM

Water-soluble and highly conductive cream for EEG-EP-EMG recordings

- Non irritating, non greasy
- Non damaging probes
- 260 g. PP-HDPE bottle (25 pcs per box)

EEG PASTE "SUPERIOR QUALITY"

Highly adhesive & above-average conductive paste, suitable for disc electrodes

- Pale white colour
- Water-soluble
- PH from 6.6 to 7
- Low impedance
- 250 g. jar (36 pieces per carton)

LUBRI SUPERGEL



LUBRI SUPERGEL

Lubricant gel for both professional and family use.

Thanks to its ingredients, this gel is the best for endoscopy, gastroscopy, colonoscopy, proctoscopy and, more in general, for all those domestic applications that do not need a sterile product.

Hypoallergenic, non-irritant, non-toxic, odourless, water soluble.

Packing:

- 82 gr. bottle ; n°48 bottles/box
- 260 gr. bottle ; n°25 bottles/box;

AQUASONIC PARKER



AQUASONIC - PARKER

Ceracarta is the official distributor of this gel for ultrasounds famous all over the world. The quality and characteristics of this gel are undoubtedly excellent.

Available in three versions:

- 250 gr (box with 12 bottles - carton with 72 pieces)
- 5 kg (individual package with bottle/refill - box with 4 pieces)
- 20 gr Sterile gel in a small metal bag wrapped inside another sealed bag - guaranteed (48 pieces per box)



CUSTOMISED GELS

CUSTOMISED GELS

Ceracarta can provide solutions for important customers who, for any reason, might require a customised gel packaging. At no extra cost you can have bottles or containers graphically designed to meet each customer's taste and requirements.



Via Secondo Casadei N.14
zona industriale Villa Selva - 47122 Forlì (Italy) tel. +39.0543.780055 - fax +39.0543.781404
www.ceracarta.it - e mail: export@ceracarta.it