



ULTRASOUND SYSTEMS & IMAGING MATERIAL

- Wide-ranging and continuous availability of original materials, both black and white and colour
- Centres for customer consulting for any clarification, information and sample supply
- Extreme flexibility in order dispatching and great rapidity in supplies
- Wide range of video printers, video recorders, monitors and video cameras of the major brands

 **CERACARTA**

...for over
50 years
CERACARTA
Italy
since...1961



SONY BLACK AND WHITE PAPERS

Four main models of rolls 110 mm. wide, between 18 and 20 meters long, from the UPP 110 S standard quality to UPP 110 HD high definition type, to the excellent UPP 110 HG up to the top UPP 110 HA version • **boxes of ten pieces and cartons of 100 pieces.**

Two models: UPP 210 SE and UPP 210 HD, measuring 210 mm x 25 m that, thanks to a plastic spacer can replace the discontinued UPP 216 SE/HD and work for loading into the UP-930/910 with no changes of settings, without forgetting the transparent film model UPT 210 BL • **boxes of five pieces and cartons of twenty pieces.**



SONY COLOUR PAPERS / FILMS & OTHER MEDIA

The complete range of original papers, in all sizes [A6/A5/A4], for all types of printers launched on the market over the past years, included the latest models UPC 24 SA/LA. Considering the large volume of sales, there is always a substantial amount of stock available which means that in a few days we can meet any request for any quantity. Furthermore, even though paper does not have an "expiry date", such rapid turnover makes it possible to always have "fresh" products of recent manufacturing.

MITSUBISHI COLOUR PAPERS



MITSUBISHI BLACK AND WHITE PAPERS

Also in the case of this important brand which has a substantial market share we have available all original paper types in black and white, from the basic KP 61B / KP 65 HM & KP 91 HG to KP 70B and KP 75HM

This range of colour papers is compact, but at the same time sufficiently complete to meet any requirement in terms of Mitsubishi printers launched on the market. Original paper in sheets or in rolls (Mitsubishi exclusive) such as CK 700/800/900 including the exceptional "four-layers" CK 700 4P up to the latest CK 30 S / L • constant and large stock availability.


VIDEOPRINTERS FOR ULTRASOUNDS



UP-897MD



UP-D711MD



UP-25MD



MITSUBISHI P.93

UP 897 MD / UPD 897 MD - The ideal choice for ultrasound systems - Black & white

UP - D711 MD - The compact, high speed printer - Small & lightweight

UP 25 MD / UPD 25 MD - Colour printer with high design - Perfect in a wide range of medical applications

P 93 / P 95 D - Black & white printers - Analogue & digital versions

SONY

Paper		printer mod
UPP 84 S	b/w mm. 84x13,5mt.	UP - D711 MD
UPP 84 HG	b/w mm. 84x12,5mt.	UP - D711 MD
UPP 110 S	b/w mm. 110x20mt.	UP 860/890/895/897
UPP 110 HD	b/w mm. 110x20mt.	UP 860/890/895/897
UPP 110 HA	b/w mm. 110x18mt.	UP 890
UPP 110 HG	b/w mm. 110x18mt.	UP 895/897
UPP 210 SE	b/w mm. 210x25mt.	UP 960/980/970/990
UPP 210 HD	b/w mm. 210x25mt.	UP 960/980/970/990
UPT 210 BL	b/w mm. 210x12,5mt.	UP 980/970/990
UPC 510	col. 200 prints A5	UP 51MD/51MDP
UPC 21 S	col. 240 prints "Small"	UP 20/21MD/D23MD/25MD/D25MD
UPC 21 L	col. 200 prints "Large"	UP 20/21MD/D23MD/25MD/D25MD
UPC 24 SA	col. 180 prints "Small"	UP 25/UPD25
UPC 24 LA	col. 160 prints "Large"	UP 25/UPD25
UPC 55	col. 200 prints A5	UP 55MD/D55
UPC-R80 MD	col. 100 prints A4	UPD - R80 MD
UPT 510 BL	Blue film 125 sheets 8"x10"	UP DF550/DF750
UPT 512 BL	Blue film 125 sheets 10"x12"	UP DF550/DF750
UPT 514 BL	Blue film 125 sheets 11"x14"	UP DF550/DF750
UPT 517 BL	Blue film 125 prints 14"x17"	UP DF500
UPT 725	b/w 100 prints 8"x10"	UP D71XR/D72XR/D74XRD
UPT 735 BL	Blue film 100 prints 8"x10"	UP D71XR/D72XR
UPT 736 BL	Blue film 100 sheets 8"x10"	UP D74XRD
UPC 740	col. 72 prints A4	UP D70A/D70APSW
UPC-770	col. 72 prints A4	UP - D 75/77 MD

MITSUBISHI

Paper		printer mod
KP 61 B e KP 61 CE	b/w mm. 110x20mt	P 60/61/65/66/67/68
		/90/91/93E/95DE
KP65HM-KP65HM-CE	b/w mm. 110x20mt	P 60/61/65/66/67/68
		/90/91/93E/95DE
KP 91 HG	b/w mm. 110x18mt	P 90/91/93E
KP 70 CE	b/w mm. 216x20mt	P 71/75/78/500E
KP 75 HM - CE	b/w mm. 216x20mt	P 71/75/78/500E
CK 700	col. mm. 110x22mt	CP700/710E
PK 700 S	col. 200 prints "Small"	CP700/710E
PK 700 L	col. 130 prints "Large"	CP700/710E
CK 700 S 4P	col. 110 prints 4 panels "Small"	CP700/710E
CK 700 L 4P	col. 75 prints 4 panels "Large"	CP700/710E
CK 900 S	col. 200 prints "Small"	CP900/910E
CK 900 L	col. 130 prints "Large"	CP900/910E
CK 30 S	col. 240 prints "Small"	CP30DW
CK 30 L	col. 200 prints "Large"	CP30DW

FUJI - ALOKA - TOSHIBA - HITACHI - KOWA - PANASONIC - ETC.

Original black and white or colour papers are available also for printers of these brands which, even though less widespread, are constantly present on the market.

ALTERNATIVE PAPERS FOR SONY - MITSUBISHI



Ceracarta uses Japanese raw materials to manufacture alternative papers in black and white of excellent quality for the models UPP 110 S, K 61B and the high definition models UPP 110 HD, UPP 110 HG, UPP 210 HD, K 65 HM, etc.

- The prices quoted are very competitive
- They are all sent in an elegant and practical package.

OTHER PAPERS AND FILMS



Brands such as **KODAK**, **CODONICS**, **AGFA**, **FUJI**, etc... are present with high-performance printers in the radiology and nuclear medicine units. Also for this important sector, Ceracarta can provide all types of papers and films.

X-RAY FILMS



A wide range of X-Ray films (mainly dry) for all the most sold brands [AGFA / FUJI / KODAK and all...]. Stock availability - Quality service.

BATTERIES



Broad range of batteries for all customers' needs both in medical and professional field. • The most popular brands along with the best quality. • Always available in stock.



Via Secondo Casadei N.14
zona industriale Villa Selva - 47122 Forlì (Italy) tel. +39.0543.780055 - fax +39.0543.781404
www.ceracarta.it - e mail: export@ceracarta.it